

1 9 in 2005 to 6 out of 6 in 2007. This  
2 includes today's and tomorrow's presentations.

3 And you will here Dr. Duggirala presenting  
4 today on behalf of the Epidemiology Branch.

5 And this is our vision. To  
6 conclude, I just would like to say that, you  
7 know, we would like to see that all important  
8 post-market questions are addressed by  
9 post-approval studies, that studies are  
10 realistic and founded on good science.

11 We would like also studies to be  
12 timely, accurate, and provide useful  
13 information, based on which we might base some  
14 of the regulatory actions if needed.

15 Also, we certainly would like to  
16 have reports that are clearly identified and  
17 effectively track. And we are committed to  
18 keep our stakeholders apprised.

19 I cannot stress enough how  
20 important it is for us to maintain and  
21 cultivate our cooperation with our pre-market  
22 colleagues because they are the technical

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114720**

A2331

1 experts for the product. And we bring the  
2 epidemiologic expertise in the study design.

3 If we proactively address all of  
4 these issues in a timely fashion during  
5 pre-market, we predict that we will have less  
6 enforcement options.

7 Just for your information, this is  
8 the Epidemiology Branch. And the current  
9 staff that is involved in cardiovascular  
10 devices are marked in blue. You will see  
11 there are three epidemiologists, one team  
12 leader, a branch chief, and the three project  
13 managers that handle post-approval study  
14 commitments with regard to cardiovascular  
15 devices.

16 Again, the post-approval studies  
17 transformation, vision, and goals present high  
18 expectations of us and of the stakeholders.  
19 Heightened expectations often bring heightened  
20 concerns about burdens, workload, perceived  
21 fairness, and added value. It is up to us and  
22 our stakeholder to discuss them openly,

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114721**

A2332

1 responsibly, and collaboratively.

2 We understand the concerns, but we  
3 have to put them into larger contexts of  
4 asking and answering the right post-market  
5 questions. We welcome an exchange of ideas on  
6 diverse methodologies that may be  
7 cost-effective, innovative, and productive.  
8 We value all analytical approaches and data  
9 sources that will give us high-quality answers  
10 to the right post-market questions.

11 Thank you.

12 CHAIRPERSON YANCY: Thank you very  
13 much. Obviously this information is important  
14 as increasingly we rely upon the repository of  
15 information from post-marketing studies to  
16 help us understand the impact of the  
17 technologies that we are considering.

18 Are there any questions for the  
19 speaker that you just heard?

20 (No response.)

21 CHAIRPERSON YANCY: Great.

22 1ST OPEN PUBLIC HEARING

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

[www.nealrgross.com](http://www.nealrgross.com)

(202) 234-4433

**CORD114722**

A2333

1 CHAIRPERSON YANCY: We will now  
2 proceed with the open public hearing portion  
3 of this meeting. Both the Food and Drug  
4 Administration and the public believe in a  
5 transparent process for information-gathering  
6 and decision-making.

7 To ensure such transparency of the  
8 open public hearing session of the Advisory  
9 Committee meeting, FDA believes that it is  
10 important to understand the context of any  
11 individual's presentation.

12 For this reason, FDA encourages  
13 you, the open public hearing, our industry  
14 speaker, at the beginning of your written and  
15 oral statements to advise the Committee of any  
16 financial relationship that you may have with  
17 the sponsor; its product; and, if known, its  
18 direct competitors. For example, this  
19 financial information may include the  
20 sponsor's payment of your travel, lodging,  
21 other expenses in connection with your  
22 attendance at the meeting.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114723**

A2334

1           Likewise, FDA encourages you at the  
2           beginning of your statement to advise the  
3           Committee if you do not have any such  
4           financial relationships. If you choose not to  
5           address this issue of financial relationships  
6           at the beginning of your statement, it will  
7           not preclude you from speaking.

8           There is no one that has signed up  
9           for this session. However, if there is anyone  
10          in the audience who would like to speak, we  
11          would appreciate hearing from you.

12                       (No response.)

13           CHAIRPERSON YANCY: Since no one is  
14           coming forward, we will proceed with today's  
15           agenda. Please note there will be a second  
16           opportunity as there is another open public  
17           session in the afternoon.

18           We will now proceed to the sponsor  
19           presentation. Whomever speaks first, if you  
20           can help me with the nomenclature, that would  
21           be great, the XIENCE V Everolimus-Eluting  
22           Coronary Stent System.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114724**

A2335

1 I would like to remind public  
2 observers at this meeting that while this  
3 meeting is open for public observation, public  
4 attendees may not participate except at the  
5 specific request of the panel.

6 We will begin with the sponsor  
7 presentation. Thank you.

8 MR. JOHNSON: Thank you, Mr.  
9 Chairman.

10 SPONSOR PRESENTATION

11 MR. JOHNSON: The proper name is  
12 the XIENCE V Everolimus-Eluting Coronary  
13 Stent.

14 CHAIRPERSON YANCY: Thank you.

15 MR. JOHNSON: Sure. Good morning.  
16 My name is Gary Johnson. I am Vice President  
17 of Regulatory Affairs, Clinical Research, and  
18 Quality Assurance for Abbott Vascular.

19 And on behalf of the employees of  
20 Abbott Vascular, our physician investigators,  
21 and the patients enrolled in our clinical  
22 studies, I would like to thank the panel

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114725**

A2336

1 members and the Food and Drug Administration  
2 for this opportunity to present the XIENCE V  
3 Everolimus-Eluting Coronary Stent program.

4 During today's presentation, we  
5 will cover a number of topics. First, I will  
6 do a brief introduction. Then Dr. Murthy  
7 Simhambhatla from Abbott will review XIENCE V  
8 design goals.

9 Dr. Leslie Coleman from Abbott will  
10 summarize our preclinical animal studies.  
11 Then Dr. Gregg Stone from Columbia University,  
12 who is the PI of the SPIRIT III clinical  
13 trial, will review the results of three  
14 randomized clinical trials: SPIRIT FIRST,  
15 SPIRIT II, and SPIRIT III.

16 Dr. Stone will be followed by Dr.  
17 Mitchell Krucoff from Duke University, who is  
18 the co-PI for XIENCE V U.S.A. post-approval  
19 clinical study. He will review a combined  
20 safety analysis and our integrated  
21 post-approval clinical strategy. Finally, Dr.  
22 Krishna Sudhir from Abbott will close and

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114726**

A2337

1 summarize our presentation today.

2 In addition to the presenters, we  
3 also have several consultants with us today:  
4 Professor Stuart Pocock from the London School  
5 of Hygiene and Tropical Medicine; Dr.  
6 Alexandra Lansky from CRF, who served as the  
7 angiographic core laboratory; Dr. Peter  
8 Fitzgerald from Stanford University, who  
9 served as the IVUS core laboratory; and Dr.  
10 Renu Virmani from CVPPath International. We  
11 also have Mr. Ron Van Valen with us from  
12 Novartis today.

13 The purpose of our presentation  
14 today is fourfold. First, we want to review  
15 XIENCE V design goals and provide a detailed  
16 understanding of the major design  
17 characteristics and why they were selected.

18 Second, we want to review the  
19 breadth and depth of our preclinical animal  
20 studies and vessel healing evaluations.  
21 Third, we want to demonstrate the XIENCE V  
22 clinical data in its totality establishes a

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114727**

A2338



1 reasonable assurance of safety and  
2 effectiveness based on valid scientific  
3 evidence.

4 And, finally, we want to review the  
5 XIENCE V post-approval clinical strategy that  
6 augments our pre-approval clinical data and is  
7 effectively powered to evaluate low-frequency  
8 events.

9 We are seeking an indication for  
10 XIENCE V which is consistent with other  
11 drug-eluting stents. The proposed indication  
12 is for improving coronary luminal diameter in  
13 patients with symptomatic heart disease due to  
14 de novo native coronary artery lesions with  
15 lengths less than or equal to 28 millimeters  
16 and with reference vessel diameters of 2.5  
17 millimeters to 4.25 millimeters.

18 We are seeking approval for five  
19 diameters of stents, 2.5, 2.75, 3.0, 3.5, and  
20 4.0 millimeters. These diameter stents will  
21 be available in six lengths: 8, 12, 15, 18,  
22 23, and 28 millimeters. These diameters and

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114728**

A2339

1 lengths will be available in both a rapid  
2 exchange and over the wire delivery system.  
3 All stem sizes have consistent drug dose  
4 density of 100 micrograms per centimeter<sup>2</sup>.

5 There are three major design  
6 components of XIENCE V coronary stent system:  
7 the stents and delivery system, which are  
8 based on the approved Multi-Link VISION and  
9 Multi-Link MINI VISION coronary stent systems;  
10 the drug matrix, which is a fluorinated  
11 copolymer that has previously been approved on  
12 other vascular application devices; and the  
13 drug, everolimus, which is manufactured by  
14 Novartis Corporation.

15 Everolimus under the brand name  
16 Certican has received two approvable letters  
17 from FDA for organ transplant indication.  
18 Novartis has also granted FDA rights to  
19 reference their IND and NDA to support XIENCE  
20 V PMA review. XIENCE V stent has received  
21 regulatory approval and is currently marketed  
22 in 64 countries outside of the United States.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114729**

A2340

1           To put our regulatory and clinical  
2           activities into perspective, I wanted to  
3           provide a brief regulatory history. Abbott  
4           Vascular worked collaboratively with FDA in  
5           late 2004 and early 2005 to develop a SPIRIT  
6           III pivotal clinical trial design.

7           At the time of initiation at the  
8           SPIRIT III clinical trial, in May of 2005, FDA  
9           agreed that a trial design and the supporting  
10          clinical data from the XIENCE V clinical  
11          program would provide adequate assurance of  
12          safety and effectiveness for the XIENCE V  
13          system.

14          During that process, FDA reviewed  
15          the everolimus safety, pharmacology,  
16          toxicology, and ADME studies and identified no  
17          concerns. FDA considers everolimus to be a  
18          well-characterized and studied drug;  
19          therefore, not a new molecular entity, or NME.

20          Since everolimus was not an NME,  
21          the requirement for 2,000 treated patients in  
22          clinical studies typically required for a

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114730**

A2341

1 drug-eluting stent with an NME did not apply  
2 to the XIENCE V program.

3 The XIENCE V's clinical program  
4 will be overviewed in detail today. It is a  
5 comprehensive, integrated pre-approval and  
6 post-approval program that includes over  
7 16,000 patients. It includes four  
8 pre-approval clinical studies, for which you  
9 will be presented today, and six ongoing or  
10 planned clinical studies, which will be  
11 reviewed in more detail later in the  
12 presentation.

13 In summary, the pre-approval  
14 clinical studies in their totality have  
15 demonstrated the following. All the trials  
16 have met their pre-specified primary and  
17 powered major secondary endpoints. They have  
18 demonstrated non-inferiority and superiority  
19 in late loss or bare metal stent. They have  
20 demonstrated non-inferiority and superiority  
21 in late loss over an approved drug-eluting  
22 stent. They have also shown non-inferiority

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114731**

A2342

1 in target vessel failure compared to the TAXUS  
2 drug-eluting stent. And, additionally, all of  
3 these trials will have long-term follow-up  
4 after five years.

5 The ongoing and planned clinical  
6 studies in their totality are designed to  
7 include real-world patients, their power to  
8 effectively detect low-frequency events of .5  
9 percent, their design to support label  
10 expansion to more complex patient subsets.  
11 And these studies will also have long-term  
12 follow-up after five years.

13 In addition to our planned  
14 analysis, in response to the panel's comments  
15 in December, we have also performed a safety  
16 subset analysis of all available two-year data  
17 from SPIRIT II and SPIRIT III.

18 Results of this analysis are  
19 consistent with the one-year data from SPIRIT  
20 II and III as well as the three-year data from  
21 SPIRIT FIRST. The results of this will be  
22 presented in more detail later in the

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114732**

A2343

1 presentation.

2 So, with that brief introduction, I  
3 would now like to turn the podium to Dr.  
4 Murthy Simhambhatla to review the XIENCE V  
5 design.

6 DR. SIMHAMBHATLA: Good morning.  
7 Work on the XIENCE V system began at a time  
8 when the first iteration of drug-eluting  
9 stents were already on the international  
10 market. It was our design objective to  
11 develop a second generation drug-eluting stent  
12 by integrating well-categorized, well-tested,  
13 and proven components into a system capable of  
14 assuring a high level of safety,  
15 effectiveness, and deliverability.

16 We made a decision early in the  
17 design process to use the Multi-Link VISION  
18 and MINI VISION systems as the platform for  
19 XIENCE V. The Multi-Link VISION and MINI  
20 VISION stents are the number one selling bare  
21 metal stents globally and in the United  
22 States. This is a proven and well-tested

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114733**

A2344

1 system.

2 The Multi-Link VISION stent is  
3 flexible and has thin struts, characteristics  
4 that we believe are particularly important in  
5 a DES era, by assuring a high level of  
6 conformability to the coronary vasculature and  
7 potentially a greater extent of endothelial  
8 cell coverage. The Multi-Link VISION stent is  
9 also proven to be highly deliverable, a  
10 characteristic that we wish to preserve, even  
11 after putting a drug coating on the stent.

12 A second design objective was to  
13 develop a thin, biocompatible drug coating.  
14 We felt that a thin coating would minimize the  
15 total cross-section of the coated stent strut  
16 and that by doing so, you could not only  
17 potentially facilitate the extent of  
18 endothelialization but also minimize the  
19 potential for flow impairment of side branches  
20 traversed by the stent struts.

21 In order to develop a thin  
22 biocompatible drug coating, the system had to

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

[www.nealrgross.com](http://www.nealrgross.com)

(202) 234-4433

**CORD114734**

A2345

1 be effective at low drug loading. The polymer  
2 had to be stable. Particularly in the local  
3 in vivo milieu, the coating had to be uniform  
4 and conformal to minimize the potential for  
5 plated adhesion to the stent surface. And the  
6 drug release had to be controlled and compete  
7 over time in order to reduce the potential for  
8 persistent vascular effects related to the  
9 drug. And, finally, the system in totality  
10 had to exhibit good hemocompatibility and  
11 vascular compatibility.

12 Shown here are the four components  
13 of the XIENCE V system, as previously  
14 described. The platform is a Multi-Link  
15 VISION stent and stent delivery system. The  
16 drug is everolimus. And the drug-carrying  
17 matrix is a fluorinated copolymer. I will now  
18 discuss each of these components in turn  
19 relative to our design objectives.

20 The Multi-Link VISION stent, which  
21 is the platform for XIENCE V, is based on  
22 cobalt chromium technology. This technology

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114735**

A2346



1 allows us to develop thinner struts without  
2 compromising radial opacity or radial  
3 strength. The thin struts also allow for good  
4 flexibility and conformability and in  
5 combination with the delivery system result in  
6 a low system profile. And, finally, the  
7 delivery system itself has been optimized and  
8 minimized vessel injury outside a stented  
9 segment by reducing the amount of balloon  
10 overhang outside the stent.

11 Shown here are in vitro flow data  
12 by Julio Palmaz and his colleagues. These  
13 data indicate that the extent of endothelial  
14 coverage is related to the barrier to flow.  
15 In particular, these data indicate that the  
16 extent of endothelial coverage is compromised  
17 for obstacle thicknesses exceeding 100  
18 microns. Based on this result, this group  
19 postulated that endothelial coverage may be  
20 imputed for thicker thin struts.

21 It is in this context that we  
22 believe that the progression toward thinner

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114736**

A2347

1 thin struts is important. Shown on the left  
2 are the scanning electron microscope  
3 cross-sections for the first iteration of  
4 drug-eluting stents, for which the strut  
5 thickness is significantly greater than 100  
6 microns. Shown on the right is the  
7 cross-section for XIENCE V stent with a strut  
8 thickness of 81 microns.

9 Also of note is the dark outline of  
10 the polymer coating around the bright stent  
11 struts. Shown on the left are the outlines  
12 for the first iteration of regulating stents  
13 where the coating thickness varies from 13 to  
14 20 microns. The coating thickness for the  
15 XIENCE V stent, on the other hand, is  
16 approximately eight microns.

17 We tested the deliverability of the  
18 XIENCE V system extensively in our synthetic  
19 coronary artery models that simulated  
20 tortuosity and simulated lesions. In these  
21 three-dimensional models, we demonstrated with  
22 the deliverability of the XIENCE V system for

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

[www.nealrgross.com](http://www.nealrgross.com)

(202) 234-4433

**CORD114737**

A2348

1 the smallest and largest stent sizes were  
2 equivalent to the MINI VISION and VISION bare  
3 metal stents, respectively.

4 The drug everolimus is a  
5 proliferation signal inhibitor that acts in  
6 the late G1 phase of the cell cycle to inhibit  
7 cellular proliferation in a reversible manner.

8 Everolimus belongs to the same family of  
9 synthetic macrolide compounds as sirolimus.  
10 And both these drugs have IC50 values in a  
11 similar range for the inhibition of smooth  
12 muscle cell proliferation.

13 We studied a wide range of drug  
14 doses in porcine coronary arteries from 100  
15 micrograms per centimeter<sup>2</sup> to 800 micrograms  
16 per centimeter<sup>2</sup>. We observed sufficient drug  
17 effect at 100 micrograms per square centimeter  
18 with no evidence of toxicity of medial  
19 necrosis at 800 micrograms per centimeter<sup>2</sup>.

20 The lowest effective dose of 100  
21 micrograms per centimeter<sup>2</sup> was there for  
22 selective clinical development. With this

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114738**

A2349

1 dose, we had an eight-fold safety margin given  
2 the absence of toxicity of medial necrosis  
3 with the 800 microgram per centimeter<sup>2</sup> dose in  
4 our porcine models.

5 The amount of drug on the XIENCE V  
6 stent is significantly reduced relative to  
7 other limus-eluting stents. In particular,  
8 the amount of drug in the XIENCE V stent is  
9 reduced by 41 percent relative to Cypher.  
10 This is notable given that the IC<sub>50</sub> values for  
11 both everolimus and serolimus are in a similar  
12 range for the inhibition of full muscle cell  
13 proliferation.

14 We, therefore, achieved a key  
15 design objective of effectiveness with reduced  
16 drug loading. And this will be demonstrated  
17 later in the clinical presentation.

18 The XIENCE V polymer selection and  
19 coating design were optimized for the  
20 controlled elution of everolimus over time and  
21 for the complete release of drug over time.  
22 Approximately 80 percent of the drug is

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114739**

A2350

1 released by 30 days. And substantially all of  
2 the drug is eluted by 820 days.

3 XIENCE V coating design comprises a  
4 primer and matrix system. In the expanded  
5 view on the right, the stent strut is shown in  
6 white. The drug-carrying fluoropolymer matrix  
7 is shown in blue. And a thin primer layer is  
8 shown in red. It is the function of the  
9 primer to ensure good adhesion between the  
10 drug coating and the thin strut.

11 This system does not have a top  
12 coat. In our experience, this system allows  
13 for better manufacturing control and drug  
14 release than a top coat system for such thin  
15 coatings. This system also allowed us to  
16 optimize the adhesion of the coating to the  
17 stent strut while minimizing unwanted  
18 adhesions to the delivery balloon.

19 The drug-carrying matrix is an  
20 ultra pure copolymer comprised of vinylidene  
21 fluoride and hexafluoropropylene monomers.  
22 This polymer has been used in approved

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114740**

A2351

1 cardiovascular, neurologic, and ophthalmic  
2 sutures.

3 The ratio of the vinylidene  
4 fluoride and hexafluoropropylene allows us to  
5 optimize the coating elasticity in order to  
6 prevent the coating from cracking upon stent  
7 expansion and coating toughness to ensure the  
8 durability of the coating during the act of  
9 stent delivery to the target lesion.

10 This polymer is one of the most  
11 stable entities chemically because of its  
12 durable carbon carbon backbone and the  
13 covalent carbon fluorene bonds. And this  
14 stability confers to this polymer a high  
15 degree of stability in vivo as well as  
16 biocompatibility. And, finally, this polymer  
17 has good hemocompatibility.

18 Shown here are micrographs of the  
19 XIENCE V system, illustrating its coating  
20 integrity. The coating was designed to  
21 minimize webbing, bridging, and strut-to-strut  
22 contact in the crimped state. It was also

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114741**

A2352

1 designed to maintain the coating integrity  
2 after simulated use, stent expansion, and  
3 fatigue testing.

4 A key design objective for the  
5 XIENCE V system was to assure a level of  
6 hemocompatibility that was at least equivalent  
7 to the bare metal VISION platform. We tested  
8 hemocompatibility in accordance to ISO10-993  
9 and showed in an un-hecronized ex vivo shunt  
10 study that the amount of polymers accumulated  
11 on the XIENCE V stent was less than that on  
12 the bare metal VISION stent. We, therefore,  
13 surpassed our objective of ensuring equivalent  
14 hemocompatibility to the bare metal stent.

15 We also studied the vascular  
16 response of the XIENCE V system and XIENCE V  
17 copolymer extensively in porcine models and  
18 demonstrated that all the way out to two  
19 years, the polymer response is equivalent to  
20 the VISION bare metal stent. We have also  
21 studied the vascular response of three times  
22 the amount of polymer on the stent and have

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114742**

A2353

1 found the response equivalent to the bare  
2 metal vision stent.

3 So, in summary, the XIENCE V system  
4 is built on the proven VISION stent and stent  
5 delivery system. The VISION stent is flexible  
6 and has thin struts. It is also a deliverable  
7 stent.

8 We have also developed a thin,  
9 biocompatible drug coating that is effective  
10 at low drug loading. The polymer is stable.  
11 The coating is uniform and conformal around  
12 the stent struts.

13 The drug release is well-controlled  
14 and complete over time. And, finally, the  
15 system exhibits good hemocompatibility and  
16 vascular compatibility.

17 I will now turn over the podium to  
18 my preclinical colleague, Dr. Leslie Coleman.

19 DR. COLEMAN: Good morning. I  
20 would like to present to you an overview of  
21 the XIENCE V preclinical program. The  
22 clinical program consisted of an extensive

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114743**

A2354



1 assessment of the biocompatibility of the  
2 XIENCE V system and assessment and  
3 characterization of the pharmacokinetics of  
4 XIENCE V, a comprehensive safety assessment,  
5 and an assessment of the endothelial cell  
6 response to XIENCE V.

7 The biocompatibility of the XIENCE  
8 V system was demonstrated through numerous in  
9 vitro and in vivo studies. All studies were  
10 conducted in compliance with applicable  
11 guidelines, and all studies passed.

12 The pharmacokinetics of the XIENCE  
13 V was characterized in a porcine coronary  
14 artery model. And, as you can see in the  
15 graph, the graph on the left demonstrates that  
16 the XIENCE V released everolimus in a  
17 consistent and controlled manner, with  
18 complete drug release by 120 days. And we  
19 believe that it's very important to have  
20 complete release of the drug from the system  
21 in order to allow for vessel healing.

22 These release kinetics translate

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114744**

A2355

1 into effective arterial delivery, as you can  
2 see on the right graph, where we have  
3 controlled release of everolimus to the target  
4 tissue or the stented artery over time. This  
5 has allowed for the presence of everolimus  
6 during the first several months following  
7 stent implantation consistent with the peak  
8 cellular phases of neointimal hyperplasia.

9 The clinical pharmacokinetics of  
10 XIENCE V were studied in several substudies  
11 within the SPIRIT II and SPIRIT III clinical  
12 trials.

13 Results from all P-K substudies  
14 were consistent across geographies and showed  
15 limited systemic exposure of everolimus up to  
16 a total dose of 588 micrograms. And at all  
17 times the amount of systemic everolimus  
18 correlated with the number of stents implanted  
19 into the patient.

20 Importantly, systemic exposure to  
21 everolimus was well below the minimal  
22 therapeutic blood level of three nanograms per

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

**CORD114745**

A2356

1 ml that must be maintained at a steady state  
2 when everolimus is delivered orally to prevent  
3 organ transplant rejection.

4 We conducted a comprehensive safety  
5 assessment of the XIENCE V system. This  
6 entailed 35 animal studies. We evaluated the  
7 system in two species. And we have data on  
8 the XIENCE V system extending from 28 days out  
9 to 2 years.

10 We chose to evaluate the XIENCE V  
11 in two animal species in order to account for  
12 any species-specific responses. And we  
13 evaluated the XIENCE V system in numerous  
14 configurations, including single stents, which  
15 we did evaluate in two species, again with  
16 data out to two years. We evaluated the  
17 response to overlapping XIENCE V stents, also  
18 in two species.

19 And then, as mentioned previously,  
20 we conducted studies to understand fully the  
21 safety margin of the system by evaluating a  
22 maximum dose system that has eight times the

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114746**

A2357

1 amount of everolimus on the system as well as  
2 polymer-only systems, ranging from 1X polymer  
3 out to three times the amount of polymer on  
4 the system. And we have evaluated these  
5 polymer systems out to two years.

6 From a safety perspective, our goal  
7 was to establish effective drug delivery with  
8 rapid vessel healing. We defined vessel  
9 healing by the following four criteria that we  
10 should have a smooth muscle cell rich  
11 neointima incorporating all stent struts as  
12 rapidly as possible. There should be minimal  
13 persistent fibrin, minimal long-term  
14 inflammation, and a rapidly endothelialized  
15 lumen.

16 These are representative histologic  
17 images of XIENCE V as compared to a VISION  
18 bare metal stent at numerous time points from  
19 28 days out to 2 years. The bar graph below  
20 demonstrates the inflammatory scores over  
21 time, again comparing XIENCE V to VISION  
22 metallic stents from 28 days out to 2 years.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114747**

A2358

1 Inflammation is scored on a scale  
2 of zero to four with a score of zero to one  
3 being considered background in this particular  
4 model, in the porcine model.

5 As you can see, the neointimal  
6 response to XIENCE was similar to a VISION  
7 metallic stent. And there is minimal  
8 long-term inflammation. From 180 days out to  
9 one year and then to 2 years, the neointimal  
10 response is stable.

11 At higher magnification, one can  
12 appreciate the cellular composition of the  
13 vessel wall in response to a XIENCE implant.  
14 The bar graph demonstrates fibrin over time.  
15 And, as we would expect, there is peri-strut  
16 fibrin at 28 days consistent with peak drug  
17 elution. But as the system elutes the drug,  
18 there is minimal to no fibrin consistent with  
19 elution of the drug and no longer drug being  
20 detected in the tissue. Again, the neointimal  
21 response is stable from six months out to two  
22 years, consistent with a healed vessel.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114748**

A2359

1                   Endothelialization in the porcine  
2                   model was complete at all time points  
3                   evaluated. To summarize the safety response,  
4                   during the active phase of drug elution, we  
5                   did observe that there was a fully  
6                   endothelialized lumen, there was neointimal  
7                   coverage of all stent struts with peri-strut  
8                   fibrin, consistent with drug elution,  
9                   inflammation comparable to VISION metallic  
10                  stent and no mineralization, no medial  
11                  necrosis, demonstrating the overall lack of  
12                  vessel toxicity.

13                  At 180 days and beyond, the phase  
14                  at which there is no longer a drug detected in  
15                  the tissue, the vessels were again fully  
16                  endothelialized and largely in a quiescent,  
17                  healed state. And so these findings we  
18                  believe are consistent with vessel healing.

19                  So, to conclude the safety  
20                  assessment, we have demonstrated safety in two  
21                  animal models with data out to two years. We  
22                  have met the goal of our DES safety program by

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114749**

A2360