

1 demonstrating these parameters, that there is  
2 a neointima incorporating all stent struts,  
3 minimal fibrin, minimal long-term  
4 inflammation, and a rapidly endothelialized  
5 lumen, consistent with vessel healing.

6 We have also been very interested  
7 in understanding the endothelial cell response  
8 to a XIENCE V system. As such, we have  
9 conducted several novel research studies in  
10 collaboration with Dr. Renu Virmani in order  
11 to understand the endothelial cell response to  
12 a XIENCE V stent in comparison to a metallic  
13 VISION stent as well as other commercially  
14 available drug-eluting stents.

15 All of these studies were conducted  
16 in a rabbit model, which allows us to really  
17 differentiate the endothelial response to  
18 various stent platforms.

19 And as part of these studies, we  
20 conducted an assessment, an in-depth  
21 assessment, of the endothelial cell coverage  
22 by qualitative and quantitative means through

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114750**

A2361

1 scanning electron micrograph evaluation. We  
2 also looked for the presence of specific  
3 endothelial cell markers through both confocal  
4 microscopy and molecular evaluation.

5 This I chose representative  
6 scanning electron micrograph of the luminal  
7 surface of stents that have been cut  
8 longitudinally following 14 days in a rapid  
9 vessel, iliac vessel.

10 And, as you can see, if we look at  
11 the XIENCE V, we have good endothelial cell  
12 coverage of the stent surface. And it's  
13 somewhat similar in coverage as the VISION  
14 metallic stent. And there is greater coverage  
15 as compared to the other commercially  
16 available drug-eluting stents.

17 We quantified the endothelial cell  
18 coverage. And the data is demonstrated here.

19 And in order to assess coverage, we looked  
20 both above the stent struts as well as in  
21 between stent struts.

22 And when we look over the stent

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114751**

A2362

1 struts, you can see that there is  
2 significantly greater coverage with the XIENCE  
3 V stent as compared to other commercially  
4 available drug-eluting stents. And the XIENCE  
5 V stent has similar coverage to a VISION  
6 metallic stent.

7 In order to understand endothelial  
8 cell integrity and functionality, we looked  
9 for the presence of two specific biomarkers.  
10 We looked for expression of platelet  
11 endothelial cell adhesion molecule, which is a  
12 membrane glycoprotein that is constitutively  
13 expressed by endothelial cells. And we also  
14 looked for production of vascular endothelial  
15 growth factor, which we believe should be  
16 down-regulated with complete  
17 endothelialization.

18 The expression of PECAM is  
19 demonstrated here. This was evaluated by  
20 immunohistochemical means. And the chemical  
21 expression is shown here. And, consistent  
22 with the endothelial data that I just showed,

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114752**

A2363

1 there is significantly greater expression of  
2 PECAM on the part of XIENCE as compared to  
3 other drug-eluting stents.

4 In terms of VEGF production, we  
5 looked for both protein levels as well as gene  
6 expression. And the data shown here are from  
7 two different sets of experiments. And the  
8 data is consistent with one another.

9 We do see that the levels of VEGF  
10 expression for XIENCE are similar to a VISION  
11 metallic stent and less than the other  
12 drug-eluting stents. And we believe this is  
13 consistent with endothelialization.

14 So, to summarize, XIENCE V  
15 demonstrated rapid re-endothelialization  
16 compared to other drug-eluting stents. And  
17 XIENCE V demonstrated enhanced endothelial  
18 cell function. And, to conclude, we do  
19 believe that rapid endothelial cell coverage  
20 and function are consistent with vessel  
21 healing.

22 I would now like to introduce Dr.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114753**

A2364

1 Gregg Stone, who will speak to you about the  
2 SPIRIT clinical program.

3 DR. STONE: Thank you, Leslie.

4 Good morning. My name is Gregg  
5 Stone. I'm an interventional cardiologist at  
6 Columbia University Medical Center. And I am  
7 here to discuss the preclinical  
8 investigational pathway that has been  
9 undertaken for the SPIRIT program of the  
10 XIENCE V stent.

11 I have also represented in the past  
12 Boston Scientific as the principal  
13 investigator of the pivotal TAXUS IV trial  
14 that led to approval of that device in the  
15 United States and have also been the principal  
16 investigator of the United States TAXUS 5 and  
17 TAXUS 5 ISR trials.

18 I currently receive research  
19 support from both Abbott Vascular and from  
20 Boston Scientific, which are the manufacturers  
21 of the two stents that I will be discussing in  
22 this presentation and which were pitted

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114754**

A2365

1 against each other in the SPIRIT II and SPIRIT  
2 III trials. And I also do work with other  
3 device companies that manufacture devices that  
4 I won't be discussing today.

5 This is the second time of four you  
6 are going to see this slide describing the  
7 16,000 patients that will be enrolled, have  
8 been rolled, and will be enrolled in the  
9 pre-approval and the ongoing and planned  
10 clinical studies for the XIENCE V stent.

11 And I am going to be specifically  
12 describing the results of four different  
13 studies. These of these were randomized  
14 trials: the SPIRIT FIRST trial, the SPIRIT  
15 II, and SPIRIT III, and then the SPIRIT III  
16 4.0-millimeter stent registry arm, which was  
17 part of the SPIRIT III randomized trial.

18 So if we first look at the SPIRIT  
19 FIRST trial, this was the first in-human use  
20 of the XIENCE V stent. This trial was  
21 performed at a time when one could still  
22 compare a drug-eluting stent versus a bare

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114755**

A2366

1 metal stent. And this compared the XIENCE V  
2 stent versus the otherwise similar bare metal  
3 VISION stent.

4 This was a relatively simple  
5 randomized controlled trial in single de novo  
6 lesions. It was performed in 60 patients in 9  
7 sites in Europe. Again, relatively short  
8 focal lesions with a reference vessel diameter  
9 of three millimeters and a lesion length of up  
10 to 12 millimeters were enrolled in this study  
11 and then randomized one to one to either the  
12 XIENCE V Everolimus-Eluting Stent versus an  
13 otherwise identical multi-length VISION bare  
14 metal stent.

15 Again, this was prospective. This  
16 was single blind and randomized. Angiographic  
17 and intravascular ultrasound was performed or  
18 intended to be performed at six months and one  
19 year in all patients with clinical follow-up  
20 performed at regular intervals up to five  
21 years.

22 I will be emphasizing the primary

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114756**

A2367

1 and the major secondary endpoints for these  
2 studies. The primary endpoint was  
3 angiographic in-stent late loss at 180 days,  
4 with a major secondary endpoint of  
5 intravascular ultrasound-determined percent  
6 volume obstruction; that is, tissue growth  
7 within the stent.

8 At 180 days, this trial, of course,  
9 was underpowered to look at clinical events.  
10 Both of these endpoints looking at the degree  
11 of neointimal growth within the stents over  
12 time were powered for superiority. That is,  
13 the XIENCE V stent had to be better than the  
14 bare metal stent. And the principal  
15 investigator was Patrick Serruys from the  
16 Thorax Center.

17 What I will be doing for all of the  
18 data that I will be showing is I will be only  
19 showing you p-values when it was either a  
20 primary endpoint or a powered secondary  
21 endpoint. Otherwise we will be displaying the  
22 results as either differences or relative

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114757**

A2368



1 risks with confidence intervals, which should  
2 be considered secondary analyses or  
3 exploratory.

4 Here you can see the primary  
5 endpoint of this trial, SPIRIT FIRST  
6 drug-eluting stent, versus bare metal stent  
7 looking at in-stent late loss. That is, from  
8 the time of the immediate post-procedure to  
9 six months later, how much tissue actually  
10 accumulated at the worst spot within the stent  
11 as determined by quantitative coronary  
12 angiography?

13 And you can see there was a marked  
14 reduction in the amount of late loss from a  
15 mean of 0.85 millimeters with the VISION bare  
16 metal stent to 0.1 millimeter with the XIENCE  
17 V stent, an 88 percent reduction, which was  
18 highly statistically significant.

19 When one looks at the major  
20 secondary endpoint of percent volume  
21 obstruction, this now looks with a more  
22 sensitive intravascular ultrasound method,

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114758**

A2369

1 looking at the percentage of tissue growth now  
2 on a volumetric basis within the stent margins  
3 that were occluding the illumene, if you will.

4 You can see that, similarly, there  
5 was a marked reduction in the percent volume  
6 obstruction from almost 30 percent of the  
7 stent being filled with tissue with the VISION  
8 bare metal stent. And this was reduced 72  
9 percent to an 8 percent volume obstruction  
10 with the XIENCE V stent, again highly  
11 statistically significant.

12 Now, I did mention that this trial  
13 was underpowered for clinical events, but, of  
14 course, these patients were followed  
15 clinically. And we have data now out to three  
16 years on the patients that were enrolled in  
17 the XIENCE V stent versus the VISION stent in  
18 SPIRIT FIRST.

19 You can see that, importantly,  
20 there were no cardiac deaths in either arm.  
21 Other events were relatively low in frequency,  
22 especially in the XIENCE V arm, in these

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

**CORD114759**

A2370

1 patients followed out to three years.

2 I will point out that there were  
3 trends towards reduced target lesion  
4 revascularization. This is the purest  
5 clinical surrogate of drug-eluting stent  
6 efficacy. This means ischemia leading to a  
7 repeat procedure due to restenosis, either at  
8 the lesion itself or at the margins out five  
9 millimeters from the lesion; the composite  
10 endpoints of major adverse cardiovascular  
11 events, which I will describe further later;  
12 and target vessel failure, also tended to be  
13 reduced with the XIENCE V stent, but, again,  
14 we weren't powered to show differences for  
15 this trial. And, perhaps most importantly,  
16 there were no cases of stent thrombosis out to  
17 three years in this small study with either  
18 the XIENCE V stent or the VISION bare metal  
19 stent.

20 So the conclusions from the SPIRIT  
21 FIRST trial were that this trial met both its  
22 pre-specified primary and major secondary

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114760**

A2371

1 endpoints, demonstrating superiority of the  
2 XIENCE V stent compared to the bare metal  
3 multi-link VISION stent in reducing late loss  
4 and percent volume obstruction.

5 We now entered a phase in clinical  
6 development where most physicians were using  
7 drug-eluting stents for the majority of  
8 patients with coronary artery disease. And it  
9 no longer became feasible to compare  
10 drug-eluting stents to bare metal stents.

11 So now we will be looking at  
12 studies comparing the XIENCE V stent to the  
13 otherwise widely utilized paclitaxel-eluting  
14 TAXUS stent. So this is DES versus DES. And  
15 the first such study, which was designed in  
16 Europe, was the SPIRIT II randomized  
17 controlled trial.

18 This was a more challenging study  
19 in which high-risk patients were enrolled.  
20 Patients were enrolled with up to two de novo  
21 lesions, rather than one, with a maximum of  
22 one lesion per epicardial vessel. And the

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114761**

A2372

1 lesions were also more challenging. They  
2 could be in smaller vessels or larger vessels,  
3 ranging from 2.5 to 4.25 millimeters and,  
4 probably even more importantly, much longer,  
5 up to 28 millimeters in length.

6 A total of 300 patients were  
7 enrolled in this trial at 28 sites outside the  
8 United States. And patients were randomized  
9 three to one to either the XIENCE V  
10 Everolimus-Eluting Stent or the TAXUS  
11 Paclitaxel-Eluting Stent. So this was also a  
12 prospective, single-blind, randomized trial.

13 Angiographic and intravascular  
14 ultrasound follow-up was performed or intended  
15 to be performed at 180 days in all patients  
16 and 2 years in approximately half the  
17 patients. Clinical follow-up was performed at  
18 regular intervals up to five years.

19 Now, the primary endpoint for this  
20 300-patient trial, in which a DES was compared  
21 to a DES, was angiographic in-stent late loss  
22 at 180 days. And this trial was powered to

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

**CORD114762**

A2373

1 demonstrate the XIENCE V compared to TAXUS was  
2 first not inferior and then if it met that  
3 endpoint to test whether or not it was  
4 superior in terms of reducing in-stent late  
5 loss.

6 The powered secondary endpoint was  
7 angiographic in-segment late loss at 180 days.

8 And this was powered for non-inferiority.  
9 And Patrick Serruys at the Thorax Center was  
10 again the principal investigator.

11 This shows now for the first  
12 primary endpoint the angiographic patient flow  
13 at 6 months of 300 randomized patients.  
14 Six-month angiographic follow-up was completed  
15 in 92 percent. And Europe is very good with  
16 angiographic follow-up.

17 And if we look at the primary  
18 endpoint of in-stent late loss, one can see  
19 that not only was the XIENCE V stent shown to  
20 be non-inferior to the TAXUS stent in terms of  
21 late loss, but it was also highly  
22 statistically superior with a marked 69

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114763**

A2374

1 percent reduction in late loss, from 0.36  
2 millimeters with TAXUS to 0.11 millimeters  
3 with XIENCE V.

4 And I will just remind you this  
5 compares to 0.8 to 1 millimeters with bare  
6 metal stents. So this expanded scale, you  
7 still see this marked reduction in late loss  
8 with XIENCE V compared to the otherwise  
9 leading TAXUS stent.

10 When one looks at in-segment late  
11 loss, this now looks not only at the biologic  
12 efficacy of the drug-eluting stent device;  
13 that is, what is going on within the margins  
14 of the stent; that is, in-stent late loss, but  
15 also takes into account the five-millimeter  
16 edges. And this can look at balloon stent  
17 mismatch issues, drug diffusion effects, and  
18 other such effects. And this is more the  
19 whole lesion itself.

20 And you can see even considering  
21 in-segment late loss, while the trial was only  
22 powered for non-inferiority, you can see there

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114764**

A2375

1 tended to be a 53 percent reduction in  
2 in-segment late loss, from 0.15 millimeters  
3 with TAXUS to 0.7 millimeters with the XIENCE  
4 V stent.

5 Looking at the intravascular  
6 ultrasound measures, again, of percent volume  
7 obstruction, this is the most commonly relied  
8 upon IVUS measure. You can see that there is  
9 a statistically significant 66 percent  
10 reduction in volume obstruction from 7.4  
11 percent with TAXUS -- and, again, this usually  
12 compares to about 30 percent with the bare  
13 metal stent -- down to 2.5 percent with the  
14 XIENCE V stent. So the DES is doing what it  
15 is supposed to be doing in terms of inhibiting  
16 tissue regrowth.

17 Now, if we look at the clinical  
18 outcomes, at one year -- and this trial,  
19 again, was not powered for clinical endpoints,  
20 but, of course, it's important to look at how  
21 the patients did -- you can see the 12-month  
22 clinical follow-up, which is the data that we

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114765**

A2376



1 currently have from this trial and has been  
2 completed in 99.3 percent of the patients.

3 First, looking at safety endpoints,  
4 stent thrombosis has been adjudicated by both  
5 the pre-specified per-protocol definition and  
6 then by now the widely used Academic Research  
7 Consortium definitions; that is, definite or  
8 probable ARC stent thrombosis. And one can  
9 see that the rates of stent thrombosis out to  
10 one year were low with the TAXUS stent --  
11 that's the red -- and also very low with the  
12 XIENCE V stent.

13 If one looks at cardiac death,  
14 there was zero percent cardiac death at one  
15 year with the XIENCE V stent versus 1.3  
16 percent with the TAXUS stent. There was a  
17 numerical trend towards less myocardial  
18 infarctions with XIENCE V compared to TAXUS.  
19 And there was a borderline statistically  
20 significant reduction in target lesion  
21 revascularization.

22 Again, this is clinical restenosis.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

**CORD114766**

A2377

1       These are real procedures due to recurrent  
2       ischemia, recurrent angina, recurrent symptoms  
3       requiring rehospitalization for repeat either  
4       angioplasty or surgery. So 6.6 percent with  
5       TAXUS and 1.8 percent with XIENCE.

6               As a result, when we start looking  
7       at composite measures of safety and efficacy  
8       -- and I am not a big fan of these composite  
9       measures because you can obviously have the  
10      components going different directions. So  
11      it's important to look at the components.

12             But the first composite measure  
13      that we often look at is MACE. And this is  
14      cardiac death, myocardial infarction, or  
15      target lesion revascularization which is  
16      relatively specific to the stent itself.

17             And here you can see that for the  
18      first time in a randomized trial, we actually  
19      saw a significant reduction in MACE with one  
20      drug-eluting stent versus another from 9.2  
21      percent, which is a relatively low number,  
22      actually, with the TAXUS stent for this type

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114767**

A2378

1 of lesion and patient mixed, but reduced to  
2 2.7 percent with the XIENCE V stent.

3 Now, the next component that we can  
4 look at of what can happen in terms of repeat  
5 vascularization is what we call target vessel  
6 revascularization remote; that is, remote from  
7 the target lesion. This could be new lesions  
8 that occur in the side branches, from guide  
9 catheter trauma, from a new distal lesion,  
10 progression of disease.

11 We wouldn't expect a drug-eluting  
12 stent to either prevent this or to improve  
13 upon it. And you can see these two  
14 drug-eluting stents had similar rates of TVR  
15 remote.

16 So if we look at the second  
17 composite measure of target vessel failure,  
18 which is somewhat more of a general measure  
19 now looking at cardiac death, MI, TLR, or TVR  
20 remote. You can see, of course, this will  
21 always dilute out a little bit of the ability  
22 to see the difference between two devices.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114768**

A2379

1 And that's what we see here.

2 Nonetheless, you could see a trend  
3 towards a 51 percent reduction in target  
4 vessel failure with XIENCE compared to TAXUS,  
5 9.2 versus 4.5 percent.

6 Thus, the conclusions from SPIRIT  
7 II are that the SPIRIT II trial met its  
8 pre-specified primary endpoint, demonstrating  
9 superiority of the XIENCE V stent compared to  
10 the TAXUS stent in reducing in-stent  
11 angiographic late loss.

12 So this brings us now to the  
13 pivotal United States-based SPIRIT III  
14 randomized controlled trial, which was  
15 designed in concert with FDA to support this  
16 pre-market approval application of the XIENCE  
17 V stent versus the TAXUS stent, again DES  
18 versus DES. And, as you'll see, this trial  
19 was designed in very similar fashion to SPIRIT  
20 II.

21 So, once again, we took patients  
22 with up to two de novo lesions with a maximum

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

**CORD114769**

A2380

1 of one lesion per epicardial vessel, with very  
2 similar reference vessel diameters and lesion  
3 lengths, as in SPIRIT II, an RVD of 2.5 to  
4 3.75 millimeters and lesion lengths up to 28  
5 millimeters. This actually matches the  
6 labeling of the TAXUS stent for use in the  
7 United States.

8 And this was a much larger study,  
9 randomizing 1,002 patients at 65 United States  
10 sites. Patients were randomized two to one to  
11 the XIENCE V stent compared to the TAXUS  
12 stent.

13 This again was a prospective  
14 single-blind, randomized trial. Angiographic  
15 and intravascular ultrasound was performed at  
16 eight months in pre-specified subsets of  
17 patients. And I'll describe this for you  
18 coming up.

19 Clinical follow-up was intended to  
20 five years at regular intervals in all  
21 patients. And that's ongoing. The primary  
22 endpoint for this trial was angiographic

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114770**

A2381

1 in-segment late loss at eight months. So we  
2 have now extended from in-stent late loss,  
3 which is easier to show because that's just  
4 looking at big differences within the stent,  
5 to now in-segment late loss, which is a more  
6 comprehensive measure taking into account what  
7 goes on at the edges, where the drugs might  
8 not be able to get to.

9 So, looking at this more  
10 comprehensive measure at eight months and this  
11 trial is powered, again, for sequential  
12 non-inferiority and superiority testing of  
13 XIENCE V versus TAXUS.

14 The first 564 patients enrolled  
15 into this trial were entered into the  
16 angiographic follow-up cohort with the  
17 patients after that enrolled in a  
18 non-angiographic follow-up cohort. There were  
19 no statistically significant differences in  
20 the baseline characteristics of the patients  
21 intended for angiographic follow-up versus  
22 those who weren't.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114771**

A2382

1 Now, the major secondary endpoint,  
2 which was actually also a co-primary endpoint  
3 because both endpoints were required to be met  
4 for regulatory approval, was the first time we  
5 have actually looked at a clinical endpoint.  
6 And that was ischemia-driven target vessel  
7 failure.

8 So this is this general measure of  
9 safety and efficacy at nine months, cardiac  
10 death, myocardial infarction, or target vessel  
11 revascularization that consists of target  
12 lesion or vascularization or TVR remote.

13 With 1,002 patients, this trial is  
14 powered for non-inferiority. And that was the  
15 regulatory burden that had to be met. The  
16 trial was not powered for superiority for  
17 target vessel failure. And it was my honor to  
18 be involved as the principal investigator of  
19 this study.

20 Now, in addition, we also wanted to  
21 look at the safety and efficacy of a  
22 4.0-millimeter stent. And we often treat

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114772**

A2383

1 large vessels with drug-eluting stents. Right  
2 now we have to take the currently available  
3 3.5-millimeter stents, put them in, expand  
4 them with a larger balloon to get them up to  
5 4.0. There is not currently an approved  
6 four-millimeter drug-eluting stent on the  
7 United States market.

8 So what was worked out in concert  
9 with FDA to evaluate this since we couldn't  
10 randomize it to another 4-millimeter  
11 drug-eluting stent was to take within these  
12 patient populations the appropriate lesions  
13 that were eligible for a 4-millimeter stent --  
14 and that's a reference vessel diameter of 3.75  
15 to 4.25 millimeters -- and to do a small  
16 registry, basically just to see if the results  
17 were consistent with the XIENCE V stent and at  
18 least not inferior to the results of the TAXUS  
19 stent in the remainder of the randomized  
20 SPIRIT III trial.

21 So this was a prospective  
22 single-blind 4.0-millimeter registry compared

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114773**

A2384



1 to the concurrent TAXUS control arm from  
2 SPIRIT III with angiographic follow-up at 8  
3 months in all patients intended with clinical  
4 follow-up ongoing to 5 years. And the primary  
5 endpoint that was agreed upon for regulatory  
6 approval was angiographic in-segment late loss  
7 at eight months powered for non-inferiority to  
8 TAXUS from SPIRIT III.

9 So if we first look at these  
10 primary angiographic endpoints, this is the  
11 randomized trial, the first 564 patients. And  
12 you can see that at 8 months, angiographic  
13 follow-up was completed in 77 percent of  
14 patients.

15 Usually in the United States, we  
16 get about 75 to 80 percent. And, in fact, we  
17 have powered this trial for 75 percent  
18 angiographic follow-up. If one extends this  
19 out to another one month, then we are up to 82  
20 percent angiographic follow-up, but this is  
21 the official formal window, with 77 percent.

22 This was the primary endpoint of

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114774**

A2385

1 SPIRIT III in-segment late loss at eight  
2 months. And, once again, one can see that not  
3 only was the XIENCE V stent compared to the  
4 TAXUS stent -- not only did it meet the  
5 primary endpoint of non-inferiority, but it  
6 was also highly statistically significant  
7 superior in terms of reducing in-segment late  
8 loss across the entire lesion and at the edges  
9 from 0.28 millimeters with TAXUS to 0.14  
10 millimeter with XIENCE V, a 50 percent  
11 reduction p-value of 0.004.

12 When one looks at the  
13 four-millimeter stent, you can see that,  
14 again, the late loss was 0.17 millimeters with  
15 this, which was shown statistically to be  
16 non-inferior to the 0.28 millimeters to the  
17 TAXUS randomized control trial. You can look,  
18 however, at the confidence intervals of the  
19 difference. And you can see that it does not  
20 overlap unity. So this actually was a  
21 reduction.

22 Now, if we go back to the

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

**CORD114775**

A2386

1 randomized control trial and first look at  
2 IVUS measures, we can see the IVUS once again  
3 supports on a volumetric basis what we saw on  
4 a spot analysis with the angiogram. And that  
5 is first looking at a slightly different  
6 endpoint of neointimal hyperplasia volume.

7 This is the amount of tissue that  
8 grows within the stent margin over the  
9 eight-month follow-up period. You could see  
10 that 21-millimeter<sup>3</sup> of tissue on average grew  
11 within the TAXUS stent versus 10.1-millimeter<sup>3</sup>  
12 in the XIENCE V stent, a statistically  
13 significant 52 percent reduction. And if one  
14 looks at, again, percent volume obstruction,  
15 you can see similar types of findings, a 38  
16 percent reduction, from 11.2 percent down to  
17 approximately 7 percent, with TAXUS versus  
18 XIENCE, respectively.

19 Now, importantly, intravascular  
20 ultrasound is a very powerful tool that allows  
21 us not to look at only the tissue inside the  
22 stent but to look at abnormal vascular

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114776**

A2387

1 responses. And what we like to see with a  
2 drug-eluting stent is that it inhibits tissue  
3 growth but doesn't do anything abnormal to the  
4 vessel wall. By angiography, we didn't see  
5 any aneurysms or ectasia in the study, but  
6 ultrasound is more sensitive than that.

7 What we can do is look at the  
8 external elastic lamina volume. This is  
9 actually the size of the entire vessel by  
10 intravascular ultrasound. And what one can  
11 see is when one looks at the XIENCE V stent  
12 patients in media after implant to eight-month  
13 follow-up by ultrasound, you can see that  
14 there is no growth in the vessel. It  
15 basically stays where you left it.

16 When one looked at the TAXUS stent,  
17 we see what we have seen in other trials. And  
18 that is an expansion or positive remodeling of  
19 the vessel. So the stent doesn't change, of  
20 course, over time, but the vessel actually is  
21 positive remodeling or expanding outward.

22 And when we want to see after this

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

[www.nealrgross.com](http://www.nealrgross.com)

(202) 234-4433

**CORD114777**

A2388

1     how this translates into the most concern that  
2     we have of incomplete stent apposition, which  
3     is acquired during follow-up -- that means a  
4     stent that was well-implanted initially at the  
5     end of the procedure but the artery, which  
6     presumably from vascular toxicity is positive  
7     remodeling, so it pulls away from the stent.  
8     This is a concern that, at least anecdotally,  
9     has been related to stent thrombosis.

10           We can see that with both stents,  
11     this was actually quite low. It occurred in  
12     1.1 percent of XIENCE V patients and 2.3  
13     percent of TAXUS patients. So what we have  
14     seen here by looking at angiography and what  
15     is supported by IVUS is a stent that leads to  
16     larger lumens compared to the other  
17     drug-eluting stent without positive  
18     remodeling; that is, without vascular  
19     toxicity, and without a risk of late acquired  
20     incomplete apposition.

21           Now let's look at the clinical  
22     follow-up in SPIRIT III of the 1,002

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

[www.nealrgross.com](http://www.nealrgross.com)

(202) 234-4433

**CORD114778**

A2389

1 randomized trials. Nine-month follow-up was  
2 completed in about 98 percent of patients and  
3 12-month follow-up in approximately 97.4  
4 percent of patients. And this was the  
5 co-primary endpoint of the trial; that is,  
6 target vessel failure at nine months.

7 Again, this was powered for  
8 non-inferiority. And the XIENCE V stent was  
9 shown to be non-inferior to the TAXUS stent  
10 for the co-primary endpoint, the first time a  
11 clinical endpoint has been pre-specified,  
12 target vessel failure from 9.7 percent with  
13 TAXUS to 7.6 percent with XIENCE V, a relative  
14 risk reduction of 21 percent, but you can see  
15 the confidence interval does cross the line of  
16 unity, so not statistically significant from  
17 superiority testing.

18 We now have follow-up to all of the  
19 patients out to one year in the SPIRIT III  
20 trial. And these are the hazard curves.  
21 Again, we have got TAXUS in red and XIENCE V  
22 in blue.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114779**

A2390

1           You can see here that there tends  
2       to be less target vessel failure in the XIENCE  
3       V arm compared to the TAXUS arm, 11.1 percent  
4       versus 8.5 percent, a relative 25 percent  
5       difference, but the p-value is .8. But this  
6       trend -- and I will get back to this -- is due  
7       to what tends to be less peri-procedural  
8       non-Q-wave myocardial infarctions very early  
9       on, with then what tends to be a little bit  
10      less ischemic target lesion revascularization  
11      later on.

12           This comes out, actually, more so  
13      when one now looks at the more stent-specific  
14      composite endpoint of major adverse  
15      cardiovascular events. Again, this is cardiac  
16      death, myocardial infarction, or target lesion  
17      revascularization that is right at the site of  
18      the stent and at the edges of the stent.

19           Here you can see the difference  
20      again in peri-procedural non-Q-wave MIs. And  
21      then you do see the curve spread over time, as  
22      I will show you later, because of less target

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

**CORD114780**

A2391